




The Go-Go Elite Traveller® Plus offers style and performance for your on-the-go lifestyle. Standard features include a delta tiller with wraparound handles, feather-touch disassembly, black,

non-scuffing tires, a front basket and an 18 AH battery pack. A 300 lb. weight capacity and a maximum speed up to 4.5 mph complement the versatility of the Go-Go Elite Traveller® Plus



Feather-Touch Disassembly













### Specifications 3W: SC53 4W: SC54

Weight Capacity	300 lbs.
Maximum Speed <sup>1</sup>	Up to 4.5 mph
Ground Clearance <sup>2</sup>	1.7" at transaxle, 0.83"at anti-tip, 2.6" at side deck
Turning Radius <sup>2</sup>	3W: 34.5" 4W: 45.75"
Overall Length <sup>2</sup>	3W: 40" 4W: 40.25"
Overall Width <sup>2</sup>	21.25"
Seat-to-Ground Height	20.5"- 22"
Seat-to-Deck Height	16.25"- 17.75"
Front Tires	3W: 2" x 8" solid 4W: 2" x 8" solid
Rear Tires	3" x 9" solid
Range Per Charge <sup>1,4</sup>	3W: Up to 14.5 miles 4W: Up to 14.5 miles
Total Weight Without Batteries	3W: 86 lbs. 4W: 87 lbs.
Heaviest Piece When Disassembled	3W: 32 lbs. (rear section) 4W: 34 lbs. (front section)
Standard Seating	Type: Foldable, compact, molded plastic Weight: 23 lbs. Material: Black Vinyl Dimensions: Width: 17" Depth: 17"
Motor Output	24-volt DC motor, sealed transaxle
Brakes	Regenerative and electromechanical
Battery Requirements <sup>3,5</sup>	Size: (2) 18 AH Weight: 14 lbs. each 29.5 lbs. per battery pack
Battery Charger	Off-board, 3.5 amp
Colors	2 sets with every scooter!

#### Blue

# HCPCS Code: K0806 (Group 2, Standard)

- 1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- 2 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- 3 AGM or gel cell type required.
- 4 Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 5 Battery weight may vary based on manufacturer.

### **Features**

- Charger port on tiller
- Convenient dual voltage charger can charge battery pack on-board or off-board
- Delta tiller with wraparound handles
- Feather-touch disassembly permits simple frame separation with only one hand and easily disassembles into five super lightweight pieces for convenient transport and storage
- Front frame-mounted seat post offers maximum stability
- Front tiller-mounted basket
- Includes two sets of easily changeable, red and blue colored panels
- Long-lasting, high-intensity, front LED headlight offers optimal pathway illumination
- Microprocessor-based controller offers optimal power management and added safety features
- Pride's exclusive black, non-scuffing tires
- Removable footboard basket for additional storage (4-wheel version only)

#### Accessories

- Additional 18AH battery box
- Cane/crutch holder (single, double)
- Cell phone holder
- Cup holder
- Forearm crutch holder
- 50" Lap belt
- Oxygen tank holder
- Rear basket
- Rearview mirror
- Saddlebag
- Safety flag
- Under seat storage (assembly required)
- Walker holder
- Weather cover
- Wishbone crutch holder
- XLR USB charger

### ption

Additional 18AH battery box

## Military Patches



 ${\sf NOTE:All\ specifications\ subject\ to\ change\ without\ notice.\ The\ information}$ contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Dimensions listed are for scooter only. Overall measurements will vary based on seating and accessory selections.







